

Application Note # CA-270358

Determination of Total FAME and Linolenic Acid Methyl Esters in Biodiesel According to EN-14103

Introduction

In order for biodiesel to be used as a motor fuel or blended with petroleum diesel, it must conform to standard specifications (ASTM D 6751 or EN-14214). There are standard GC methods in use today to determine if biodiesel conforms to the standard specifications, one of which is EN-14103, used to determine the ester and linoleic acid methyl ester content. Other methods include EN-14105 / ASTM D 6584 (free and total glycerine and mono, di and triglyceride content) and EN-14110 (residual methanol). Bruker has designed GC solutions for each of these standard methods. This paper describes work conducted on the Bruker 450 gas chromatograph to analyze biodiesel (B-100) in accordance with EN-14103. EN-14103 is used to verify that the ester content of Fatty Acid Methyl Esters (FAME) is greater than 90 % (m/m) and that the linolenic acid content is between 1 % (m/m) and 15 % (m/m) consistent with the EN-14214 specifications. This method is suitable for FAME which contains methyl esters between C14 and C24.

Instrumentation

Bruker 450-GC

- Injector: Split / splitless 1177, full EFC control
- Detector: FID, full EFC control
- Automatic Liquid Sampler; Model 8410 Software
- GC control and data handling: Bruker Galaxie™ Software

Materials and Reagents

Column

- Bruker Select™ Biodiesel for FAME, 30 m x 0.32 mm x 0.25 µm, part no. CP-9080

Sample Preparation

Accurately weigh approximately 250 mg of sample in a 10 mL vial, then add 5 mL of methyl heptadecanoate solution (10 mg / mL) using a pipette.

Injection volume: 1 µL

Conditions

Injector: 250 °C, split 100 mL / min

Detector: 250 °C, FID

Carrier gas: Helium, 12 psi (83 kPa)

Oven: 210 °C isothermal




Figure 1: Biodiesel chromatogram.

Peak Identification:

C16:0	Palmitic acid
C16:1	Palmitoleic acid
C17:0	Heptadecanoic acid (int. st.)
C18:0	Stearic acid
C18:1	Oleic acid
C18:2	Linoleic acid
C18:3	Linolenic acid
C20:0	Arachidic acid
C20:1	Gadoleic acid
C22:0	Behenic acid
C22:1	Erucic acid
C24:0	Lignoceric acid
C24:1	Nervonic acid

Results and Discussion

When the conditions outlined above are applied, a chromatogram is obtained, as in Figure 1.

Calculations

The ester content (C), expressed as a mass fraction in percent, is calculated using the following formula:

$$C = \frac{(\Sigma A) - AEI}{AEI} \times \frac{CEI \times VEI}{m} \times 100 \%$$

ΣA is the total peak area from the FAME C14:0 to C24:1

AEI is the peak area of methylheptadecanoate

CEI is the concentration, in mg / mL, of the methylheptadecanoate solution

VEI is the volume, in mL, of the methylheptadecanoate solution

m is the mass, in mgr, of the sample

The linolenic acid methyl ester content (L), expressed as a mass fraction in percent, is calculated using the following formula:

$$L = \frac{AL}{(\Sigma A) - AEI} \times 100 \%$$

ΣA is the total peak area from the FAME C14:0 to C24:1

AEI is the peak area of methylheptadecanoate

AL is the peak area of linolenic acid methyl ester

Table 1: Analysis results of biodiesel.

	Area ($\mu V \cdot \text{Min}$)	Quantity (mass %)
FAME content	103139	96.6
Linolenic acid	7599.2	7.1

Table 2: Repeatability results.

	FAME (mass %)	Linolenic acid (mass %)
Average	96.4	7.1
Stand. Dev.	0.20	0.015
Rel St.Dev. (%)	0.21	0.21

Results

The biodiesel sample tested was shown to be in accordance with the requirements stated in EN-14214 (the method requirements: FAME content > 96.5 % (m/m) and linolenic acid content < 12 % (m/m)). In order to verify the integrity of the system, repeatability was determined. One sample is analyzed 15 times. See Table 2 and Figure 2. A relative standard deviation of 0.21 % was achieved. Figure 2 shows the mass % results of the subsequent injections and the absolute difference obtained compared to specification limits. In the method, the absolute difference between two test results must be:

FAME > 1.6 % (m/m)

Linolenic acid > 0.1 % (m/m)

All results obtained are within the limits specified in the method.

Conclusion

The data presented here clearly shows the applicability of the Bruker 450-GC system for the analysis of biodiesel as per the EN-14103 standard method, and does so in a repeatable manner. The biodiesel sample tested in this application note fulfilled the requirements stated in EN-14214.




Figure 2: Repeatability figures. Red lines indicate the maximum and minimum allowed variation limits specified in the method.

References

EN-14103. Fat and oil derivatives – Fatty Acid Methyl Esters (FAME) – Determination of ester and linolenic acid methyl ester contents.

EN-14214:2003. Automotive fuels – Fatty Acid Methyl Esters (FAME) for diesel engines – requirements and test methods.

Keywords
Biodiesel
EN-14103
EN-14214
FAME

Instrumentation & Software
Bruker 450 Gas Chromatograph
Bruker Galaxie™ Software

For research use only. Not for use in diagnostic procedures.